

M259B SWEATING GUARDED HOTPLATE



The following tests can be carried out

- Thermal resistance.
- Water- vapour resistance
- Water -vapour permeability
- Water-vapour permeability index

International standards applied

Thermal resistance	→	ISO 11092:1993(E) BS EN 31092:1999 IS: 2702-1965
Water-vapour resistance	→	ISO 11092:1993(E) BS EN 31092:1994
Thermal and Evaporative resistance	→	ASTM F 1868-02

Materials that can be tested

- Medical textiles product such as wound dressings, surgical gown
- Fabrics for use in clothing
- Films and Coatings
- Foams and leather(including multi layer assemblies)
- Quilts and sleeping bags
- Upholstery and similar textile or textile-like products